

BD PureWick™ System

Fax the completed form to 800-376-4913

Patient

Name: _____

DOB: _____

Primary Phone: _____ Type: _____ Contact: _____

Address: _____ City: _____ State: _____ Zip: _____

PureWick™ System

The PureWick™ System works outside the body to draw voided urine away, helping keep skin dry.



PureWick™ System for women

(includes four components: The PureWick™ Urine Collection System, PureWick™ Flex Female External Catheter, PureWick™ Urine Collection Canister and Tubing).

PureWick™ System for men

(includes four components: The PureWick™ Urine Collection System, PureWick™ Male External Catheter, PureWick™ Urine Collection Canister and Tubing).

My patient is interested in receiving information about the PureWick™ System.

Clinician

Clinician: _____ Phone: _____ Fax: _____

NPI: _____ Credentials: _____ Document Generated By: _____

Facility name: _____ Date: _____

Address: _____

City: _____ State: _____ Zip: _____

Indication: The PureWick™ Flex Female External Catheter is intended for non-invasive urine output management in adult users with female anatomy, for conditions such as urinary incontinence.

Contraindication: Do not use the PureWick™ Urine Collection System with PureWick™ Flex Female External Catheters on individuals with urinary retention.

Warnings: The PureWick™ Flex Female External Catheter is designed for single-use only. Trying to wash, sanitize or reuse it may lead to risk of infection or illness. This may damage the product and materials. The product may not work and cause injury or illness. To avoid potential skin injury, never push or rub the product against the skin during placement or removal. Do not block airflow to the product (e.g., blocking the vent hole, bedpan, barrier creams). Do not insert the product into vagina, anus, or other body orifice. Stop use if an allergic reaction occurs. After use, this product may be a potential biohazard. Dispose of in accordance with applicable local, state, and federal laws and regulations. Do not use on users with frequent episodes of stool incontinence without a fecal management system in place. Do not use on users with skin irritation or breakdown at the site. Use with caution on users who had recent surgery of the external female anatomy. Do not use on users with moderate/heavy menstruation who cannot use a tampon.

Indication: The PureWick™ Male External Catheter used with the PureWick™ Urine Collection System is intended for non-invasive urine output management in male patients.

Contraindication: Do not use the PureWick™ Urine Collection System with PureWick™ Male External Catheters on individuals with urinary retention.

Warnings: Do not use on irritated skin. Do not cover fresh surgical wounds with the device. Do not use the BD PureWick™ Male External Catheter with any material or in any position that blocks the vents on the front of the device or does not allow for airflow through the device. To avoid potential skin injury, never pull the BD PureWick™ Male External Catheter direct away from the user. Always peel in the direction from head to foot. Stop use if an allergic reaction occurs. After use, this product may be a potential biohazard. Dispose of in accordance with applicable local, state and federal laws and regulations.

Please consult package insert for more detailed safety information and instructions for use.

